

1 EEG structure

The data is organized in the EEGLab format. The different fields of the structure EEG are:

EEG.setname

EEG.filename

EEG.filepath

EEG.subject

EEG.group

EEG.condition

EEG.session

EEG.comments

EEG.nbchan Number of channels

EEG.trials Number of trials (epochs)

EEG.pnts Number of samples

EEG.srate Sampling rate (Hz)

EEG.xmin Minimum time

EEG.xmax Maximum time

EEG.times Time. It is expressed in samples for non epoched data and in milliseconds for epoched data. Size = 1 x pnts

EEG.data Matrix with the data. Size = nbchan x pnts x trials

EEG.icaact

EEG.icawinv

EEG.icasphere

EEG.icaweights

EEG.icachansind

EEG.chanlocs Localization of the channels

EEG.urchanlocs Original channel localization

EEG.chaninfo

EEG.ref

EEG.event Structure with the events in the present data. It has fields

`.event(i).type` Name of each event

`.event(i).latency` Latency of each event (samples). If the data is not epoched it is expressed in samples in terms of the whole data. If the data is epoched the samples are counted as if the data were reshaped by concatenating the epochs along the time dimension.

`.event(i).urevent` Original event number. Index to recover original information regarding the event in EEG.urevent.

`.event(i).info` Other fields can be present if extra information was recovered from the event file

EEG.urevent Original events. Structure with all the original events. It has fields

`.urevent(i).type` Name of each event

`.urevent(i).latency` Latency of each event (samples). If the data is not epoched it is expressed in samples in terms of the whole data. If the data is epoched the samples are counted as if the data were reshaped by concatenating the epochs along the time dimension.

EEG.eventdescription

EEG.epoch Cell array, one cell per epoch. Each cell contains the events during that epoch. It has fields

`.epoch{i}.eventtype` Name of each event

`.epoch{i}.eventlatency` Latency of each event (samples). If the data is not epoched it is expressed in samples in terms of the whole data. If the data is epoched the samples are counted as if the data were reshaped by concatenating the epochs along the time dimension.

`.epoch{i}.eventurevent` Original event number. Index to recover original information regarding the event in EEG.urevent.

`.epoch{i}.eventinfo` Other fields can be present if extra information was recovered from the event file

EEG.epochdescription

EEG.reject rejection structure created by eeglab

EEG.stats

EEG.specdata

EEG.specicaact

EEG.splinefile

EEG.icasplinefile

EEG.dipfit

EEG.history

EEG.saved

EEG.etc.eeglabvers

EEG.datfile

EEG.eventFile

EEG.artifacts structure containing different fields regarding the artifacts identified in the data

EEG.reference data used as reference. Size 1 x pnts x trials

EEG.description structure containing the mean and standard deviation for each channel

2 Import the data

- **eega_importraw** Imports the .raw data to the EEGLab format. If the file with the channels localization is provided this information is added to the EEG structure.

Import Data 1: *eega_importraw*

```
EEG = eega_importraw(filename, filechanloc)
```

Inputs:

`filename` string with the full path and file name.
`filechanloc` string with the full path and file of the channels localization

Outputs:

EEG EEGlab structure.

- **eega_addchinfo** The localization of the channels can also be added at any other moment of the analysis by using this function.

Import Data 2: *eega_addchinfo*

```
EEG = eega_addchinfo(EEG, filechanloc)
```

Inputs:

EEG EEGlab structure
`filechanloc` string with the full path and file of the channels localization

Outputs:

EEG EEGlab structure.

- **eega_infoevents** Extra information of the events can be imported from the event file at any moment using this function.

Import Data 3: *eega_infoevents*

```
EEG = eega_infoevents( EEG, pathIn, eventtype )
```

Inputs:

EEG EEGlab structure
`pathIn` string with the full path and file to the events file (.txt)
`eventtype` cell array with the names of the events for which extra information will be recovered

Outputs:

EEG EEGlab structure.

- **eega_infomatfile** Extra information can be imported from a matlab file at any moment using this function.

Import Data 4: *eega_infomatfile*

```
EEG = eega_infomatfile( EEG, pathIn, varnames )
```

Inputs:

EEG EEGlab structure
`pathIn` string with the full path and file to the events file (.mat)
`varnames` cell array with the names of the variables to recover

Outputs:

EEG EEGlab structure.

3 Check the events

4 Filtering

High pass filtering or a notch filter have to be applied in order to remove line noise (50Hz). A low pass filter is necessary to remove slow drifts in the signal. Usually data is filter above 0.2-0.5Hz. However, with data from neonates I have observed that filtering at 0.5 kills part of the signal.

- **eega_filter** Filters the data and if required subtracts the mean per electrode.

Import Data 5: *eega_filter*

```
EEG = eega_filter( EEG, 'zeromean', zeromean, 'fcutoff', fcutoff, 'ftype', ftype, 'wtype', wtype, 'filtorder', filtorder, 'minphase', minphase )
```

Inputs:

EEG EEGlab structure
zeromean (Optional - default 0)
fcutoff Frequency limits (Hz) for the filter. One or two numbers. (Optional - default [0.2 40])
ftype (Optional - default 'bandpass')
wtype (Optional - default 'blackman')
filtorder (Optional - default $\text{ceil}(3.3 / (0.2 / 250)) / 2 * 2$)
minphase (Optional - default 'false')

Outputs:

EEG EEGlab structure.

5 Motion artifacts detection

Different algorithm are applied in order to identify artifacts. The data identified as containing artifacts is stored in EEG.artifacts, which has fields:

EEG.artifacts.BCT Logical matrix of size NumChan x NumPnts x NumTrials, indicating the data detected as bad.

EEG.artifacts.BC Logical matrix of size NumChan x 1 x NumTrials, indicating the channels detected as bad per epoch.

EEG.artifacts.BT Logical matrix of size 1 x NumPnts x NumTrials, indicating the time points detected as bad per epoch.

EEG.artifacts.CCT Logical matrix of size NumChan x NumPnts x NumTrials, indicating the data interpolated.

EEG.artifacts.summary Table with the summary of the amount of data rejected (BCT), channels rejected (BC), time points rejected (BT), and data corrected (CCT).

EEG.artifacts.algorithm.stepname Names of each of the algorithms applied.

EEG.artifacts.algorithm.parameters Parameters used for each of the algorithms applied.

EEG.artifacts.algorithm.rejxstep Amount of data rejected in each rejection step.

Each algorithm gives as output a variable BCT containing the data rejected by that algorithm, and updates BCT in EEG.artifacts, such as it contains the old and new rejected data. Each algorithm rejects data when a certain measure is above or/and below a threshold. The threshold for rejection can be set in three different ways.

1. An absolute threshold is provided.
2. The threshold is defined based on the distribution of specific measure along the whole recording for each electrode.
3. The threshold is defined based on the distribution of specific measure along the whole recording for all the electrodes

Thresholds defined relative to the distribution are estimated as $\text{thres} = M \pm n \cdot \text{IQ}$, where M is the median of the distribution and IQ is the interquartile range (percentile 75 - percentile 25). If required data can also be z-score before the detection is performed. In that case the mean and standar deviation for each electrode will be used. Notice that the data has to be filtered (low pass filter) in order for the algorithms to work properly

The different motion artifact interpolation algorithms are run by the function *eega.tArtifacts*. The parameters to run the algorithms are provided to the function as an input in a structure. An example of a structure with the parameters can be created by using the functions *Example_Art1* and *Example_Art2*. The *Art* structure specifies:

Art.KeepRejPre Defines what it is done with previously rejected data.

- 1 Keep data previously rejected as rejected. Default.
- 0 Set previously rejected data as good.

Art.KeepRejCause Keep or not the information regarding which data was rejected by each algorithm.

- 1 Save which data samples were rejected by each algorithm.
- 0 Do not save. Default

Art.Filter.do Specify if the data should be filtered before detecting artifacts.

- 1 Filter.
- 0 Do not filter. Default

Art.Filter.frq If Art.Filter.do = 1 it species the frequency band to filter the data in Hz. Default [0.2 40].

Art.Detect.MaxLoop Maximum number of rejection loops.

Art.Detect.Tolerance Maximum tolerance for the new data rejected (%).

Art.Detect.Ref.loop Indexes indicating for which loops the data has to be referenced to the mean before running the algorithms. Default = [].

Art.Detect.Ref.baddata Indicates what to do with the bad data when calculating the reference. It can be: 'none' / 'replacebynan' / 'replacebymeans'. Default = 'none'.

Art.Detect.steps{i} Name of the algorithm to apply. The algorithm to apply independently are:

- *eega.tRejSpectrum*
- *eega.tRejTimeVar*
- *eega.tRejAmp*
- *eega.tRejAmpDiff*
- *eega.tRejFastChange*

- *eega_tRejNetStation*
- *eega_tRejAmpElecVar*
- *eega_tRejAmpEpochVar*

The algorithms which reject data based on the data already rejected by the previous algorithms are:

- *eega_tMask*
- *eega_tRejShortGood*
- *eega_tRejShortBad*
- *eega_tRejSmplPercCh*
- *eega_tRejChPercSmpl*

Art.Detect.loop{i} Vector with indexes indicating in which loops the algorithm i should be applied.

Art.Detect.P{i} Parameters for the algorithm i.

Depending of the algorithm it can contain different fields. Some fields are the following:

- .dozscore** Z-score or not the data.
- .twdur** For algorithms which are applied on a time window, duration of the time window (s).
- .twstep** For algorithms which are applied on a time window, step for the sliding time window (s).
- .thresh** Values for the threshold. If absolute, absolutes values, if relative number of IQs from the median.
- .relative** Set relative thresholds (1) or absolute thresholds (0).
- .selectrode** If relative thresholds, estimate it per electrode or for all the electrodes.

Bellow all the functions and their inputs and outputs are described.

- **eega_tArtifacts** This function performs the rejection based on the parameters provided in the structure *Art*. The function does loops of rejection. In each loop, it first implements independently a series of algorithm to detect artifacts. Afterwards, another set functions can be applied in order to add or remove some rejected data based on the data already data. These second set of functions are applied in the order specified in *Art*. After each loop the amount of new data rejected is computed. Rejection loops keep going till when a maximum number of loops is reached or till when the amount of rejected data is under a tolerance threshold. Before starting the rejection the data can be filtered (this is useful for cases in which the original data is not filter and not filtered data is required). For each rejection loop it is possible to specify if the rejection will be done on the original data or on the data referenced to the mean.

NOTE: the data saved in EEG.data will alway be the original data, NO filtered and NO reference.

Artifacts Detection 1: *eega_tArtifacts*

```
[ EEG ] = eega_tArtifacts( EEG, Art )
```

Input

EEG EEGlab structure.

Art structure with the parameters to perform the rejection

Output:

EEG EEGlab structure with EEG.artifacts updated.

- **eega_tRejSpectrum** Reject epochs by electrodes when the power a given frequency band is too high or too low relative to the total power. The ratio is log transformed (decibels).

Artifacts Detection 2: *eega_tRejSpectrum*

```
[ EEG, BCT, T ] = eega_tRejSpectrum( EEG, [-4 4], 'dozscore', 0, finterval, [20 40],  
    'relative', 1 )
```

Inputs:

EEG EEGlab structure

thresh Threshold for each frequency band. Size n x 2. The first column is the threshold for the lower limit and the second column for the upper limit.

dozscore Z-score or not the data. (Optional - default 0).

finterval Frequency bands to check. Matrix of size n x 2, where n is the number of frequency bands. (Optional - default [20 40]).

relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

T thresholds for each electrode and each frequency band.

- **eega_tRejTimeVar** Rejects samples in a time window when the variance across samples is bigger or lower than a threshold

Artifacts Detection 3: *eega_tRejTimeVar*

```
[ EEG, BCT, T ] = eega_tRejTimeVar( EEG, [-Inf 4], 'dozscore', 0, 'twdur', 0.5,  
    'twstep', 0.25, 'relative', 1, 'xelectrode', 0 );
```

Inputs:

EEG EEGlab structure

thresh threshold for the variance across samples. The first column is the threshold for the lower limit and the second column for the upper limit.

dozscore Z-score or not the data. (Optional - default 0).

twdur time window duration (s). (Optional - default 0.5).

twstep time window step (s). (Optional - default 0.25).

relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).

xelectrode Estimate thresholds independently for each electrode (1) or an unique threshold for all the electrodes (0).(Optional - default 1).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

T thresholds for each electrode.

- **eega_tRejAmp** Rejects samples with a too high (or low) amplitude.

Artifacts Detection 4: *eega_tRejAmp*

```
[ EEG, BCT, T ] = eega_tRejAmp( EEGout, [-4 4], 'dozscore', 0, 'relative', 1,
'xelectrode', 0 );
```

Inputs:

EEG EEGlab structure
 thresh threshold for the amplitude. The first column is the threshold for the lower limit and the second column for the upper limit.
 dozscore Z-score or not the data. (Optional - default 0).
 relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).
 xelectrode Estimate thresholds independently for each electrode (1) or an unique threshold for all the electrodes (0).(Optional - default 1).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.
 BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.
 T thresholds for each electrode.

- **eega_tRejAmpDiff** Rejects samples when the change in amplitude is too high or low (first derivative).

NOTE: This algorithm is working fine for identifying very sharp changes in amplitude, like the heart beats or jumps in the signal. This type of artifacts are successfully corrected by the target PCA algorithm. See below.

Artifacts Detection 5: *eega_tRejAmpDiff*

```
[ EEG, BCT, T ] = eega_tRejAmpDiff( EEGout, [-4 4], 'dozscore', 0, 'relative', 1,
'xelectrode', 0 );
```

Inputs:

EEG EEGlab structure
 thresh threshold for the change in amplitude. The first column is the threshold for the lower limit and the second column for the upper limit.
 dozscore Z-score or not the data. (Optional - default 0).
 relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).
 xelectrode Estimate thresholds independently for each electrode (1) or an unique threshold for all the electrodes (0).(Optional - default 1).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.
 BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.
 T thresholds for each electrode.

- **eega_tRejFastChange** Rejects samples in a brief time window when the maximum change is above a threshold. The time window is moved sample by sample.

NOTE: This algorithm is slow, so it may be better to avoid.

Artifacts Detection 6: *eega_tRejFastChange*

```
[ EEG, BCTS, T ] = eega_tRejFastChange( EEG, 4, 'dozscore', 0, 'twdur', 0.025, 'tmask', 0, 'relative', 1 )
```

Inputs:

EEG EEGlab structure
thresh threshold for the variance across samples. The first column is the threshold for the lower limit and the second column for the upper limit.
dozscore Z-score or not the data. (Optional - default 0).
twdur time window duration (s). (Optional - default 0.025).
tmask length of a mask to reject samples around the samples detected as bad (s). (Optional - default 0).
relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.
BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.
T thresholds for each electrode.

- **eega_tRejNetStation** Implements the NetStation algorithm in a sliding time window. The algorithm calculates a fast running average and a slow running average in the time window, and reject samples in the time window if the fast running average is above the threshold, or if the difference between the fast and slow running average is above the threshold.

NOTE: Do not apply this algorithm with relative thresholds per electrode because the data from some electrodes is rejected even if there are not artifacts. This was happening even if the data was referenced to the mean.

Artifacts Detection 7: *eega_tRejNetStation*

```
[ EEG, BCT, BCT_f, BCT_d, T ] = eega_tRejNetStation( EEG, 4, 4, 'dozscore', 0, 'twdur', 0.5, 'twstep', 0.25, 'relative', 1, 'selectrode', 0 )
```

Inputs:

EEG EEGlab structure
thresh_fa threshold for the fast algorithm.
thresh_da threshold for the difference algorithm.
dozscore Z-score or not the data. (Optional - default 0).
twdur time window duration (s). (Optional - default 0.5).
twstep time window step (s). (Optional - default 0.25).
relative Calculate relative (1) or absolute (0) thresholds. (Optional - default 0).
selectrode Estimate thresholds independently for each electrode (1) or an unique threshold for all the electrodes (0).(Optional - default 1).

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

- BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.
- BCT_f logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the fast algorithm.
- BCT_s logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the difference algorithm.
- T thresholds for each electrode for the fast and difference running average.

- **eega_tRejAmpElecVar** Rejects samples when the amplitude is too high relative to the amplitude in the other electrodes. $M \pm n \cdot IQ$, where the median and IQ are calculated across electrodes for each time point.

Artifacts Detection 8: *eega_tRejAmpElecVar*

```
[ EEG, bcts ] = eega_tRejAmpElecVar( EEGout, [-4 4] );
```

Inputs:

- EEG EEGlab structure
- thresh threshold for the amplitude in number of interquartile ranges.

Outputs:

- EEG EEGlab structure with EEG.artifacts updated.
- BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

- **eega_tRejAmpEpochVar** Rejects samples when the amplitude is too high relative to the amplitude in the other epochs. $M \pm n \cdot IQ$, where the median and IQ are calculated across epoch for each time point and electrode.

Artifacts Detection 9: *eega_tRejAmpEpochVar*

```
[ EEG, bcts ] = eega_tRejAmpEpochVar( EEGout, [4 4] );
```

Inputs:

- EEG EEGlab structure
- thresh threshold for the amplitude in number of interquartile ranges.

Outputs:

- EEG EEGlab structure with EEG.artifacts updated.
- BCT Logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

- **eega_tMask** Marks as bad samples around a bad segment.

NOTE: The masking instead of being done during the interpolation can also be done when artifacts are corrected.

Artifacts Detection 10: *eega_tMask*

```
[ EEG, BCTS ] = eega_tMask( EEG, 0.015 );
```

Inputs:

EEG structure.

tmask time in seconds to reject before and after a bad segment.

Outputs:

EEG EEGLab structure with EEG.artifacts updated.

BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

- **eega_tRejShortGood** Marks as bad short segments between two bad segments.

Artifacts Detection 11: *eega_tRejShortGood*

```
[ EEG, BCT ] = eega_tRejShortGood( EEG, 0.05 );
```

Inputs:

EEG structure.

tmin minimum length in seconds required of a good segment to be kept as good (s).

Outputs:

EEG EEGLab structure with EEG.artifacts updated.

BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

- **eega_tRejShortBad** Removes from the rejection matrix, BCT, segments that are too short.

NOTE: If very short artifacts (like jumps in the signal or heart beat) want to be detected, then it is better no to use this and accept even very brief artifacts.

Artifacts Detection 12: *eega_tIncShortBad*

```
[ EEG, RCT, R ] = eega_tRejShortBad( EEG, 0.05 );
```

Inputs:

EEG EEGLab structure.

minlength minimum length in seconds required of a bad segment to be kept as bad.

Outputs:

EEG EEGLab structure with EEG.artifacts updated.

RCT logical matrix of size nbchan x pnts x trials, indicating the data re-included by the algorithm.

R amount of data re-included by the algorithm.

- **eega_tRejSmplPercCh** Rejects all the electrodes at a given time point if too many channels were rejected.

Artifacts Detection 13: *eega_tRejSmplPercCh*

```
[ EEG, BCT, isbadsmpl ] = eega_tRejSmplPercCh( EEG, 0.5 );
```

Inputs:

EEG EEGLab structure.

maxBadCh maximum proportion of bad channels.

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

isbadsmpl logical matrix of size 1 x pnts x trials, indicating the time points detected as bad by the algorithm.

- **eega_tRejChPercSmpl** Rejects a channel in an epoch if too many samples were rejected.

Artifacts Detection 14: *eega_tRejChPercSmpl*

```
[ EEG, BCT, isbadch ] = eega_tRejChPercSmpl( EEG, 0.25 );
```

Inputs:

EEG EEGlab structure.

maxBadSmpl maximum proportion of bad samples.

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

BCT logical matrix of size nbchan x pnts x trials, indicating the data detected as bad by the algorithm.

isbadch logical matrix of size nbchan x 1 x trials, indicating the channels detected as bad by the algorithm.

6 Artifacts Correction and interpolation

The corrections and interpolation involves

1. To define channels that are not working properly during a whole epoch (or during all the recording time). This is done based on the amount of artifacts identified for each channel. It means defining BC.
2. To define periods with strong artifacts. This is done based on the amount of channels containing artifacts and the duration of the artifacts. It means defining BT.
3. To interpolate / correct the data when it is possible. The corrected

The correction is done in an iterative way. When it is possible, detected artifacts are corrected, and afterwards the detection is run again. The process can be repeated several times.

The correction and interpolation of artifacts can be done in two basic ways:

Target PCA During periods containing artifacts we can assume most of the variance comes from the artifacts. Under this assumption we can estimate principal components associated with artifacts during this periods, and this components can be subtracted to the real data in order to correct the artifacts. In practice, bad data segments are concatenated, and PCA is applied only on this subset of data. Then, the components explaining most of the variance are subtracted to the real data subset. Finally, corrected segments are splice into the whole data and aligned to it to avoid discontinuities. It has to be noticed that the variance of electrodes it is normalized before the method is applied, in order to have equal variability across them. Data is re-scaled at the end. Moreover, when the bad segments are concatenated they are previously de-mean.

This methods is applicable for short segments of data (no much longer than 100ms) affected by artifacts, an in particular when they have a specific spatial structure. It is particularly good for removing heart beats and jumps in the signal.

Spatial Interpolation Data from some electrodes is interpolated using the data from other electrodes by using a spherical spline. It is applicable when a few electrodes are not working properly.

For both methods it is necessary to low pass filter the data after they are applied (on non epoched data), because drifts in the signal appear after splicing the corrected data segments.

And the set of functions to do so are:

- **eega_tTargetPCAEEG** This function corrects bad segments of data using target PCA. PCA is applied only on the bad. The selection of the periods to apply PCA is made based on the length of the artifacts and the amount of bad electrodes.

NOTE: The selection of the bad segments is not trivial, the next function deals with this issue by applying the method recursively by electrode.

Artifacts Correction 1: *eega_tTargetPCAEEG*

```
EEG ] = eega_tTargetPCAEEG( EEG, nSV, vSV 'maxtime', maxtime, 'mask', mask,  
'alltime', alltime, 'allchannels', allchannels, 'savecorrected', 1 );
```

Inputs:

- EEG EEGlab structure.
- nSV Number of principal components to remove from the data.
- vSV Amount variance to remove from the data (it removes the first n components of the data that removes a fraction of the variance up to vSV).
- maxtime Maximum duration of a segment in order to apply the correction (in seconds). Default 0.200.
- mask The correction is also applied on segments around bad periods (in seconds). Default 0.
- alltime Logical value defining if the correction is applied on all samples which are bad in BCT (1) or only in that ones defined as bad in BCT but for which the an artifact was defined in BT (0). Default 0.
- allchannels Logical value defining if the correction is applied on all channels (1) or only in channels defined as good in BC (0). Default 0.
- savecorrected Logical value defining if the corrected data will be marked as good in BCT (1) or not (0). Default 1.

Outputs:

- EEG EEGlab structure with the update rejection structure.

- **eega_tTargetPCAxEIEEG** This function corrects bad segments of data using target PCA. PCA is applied only on the bad data and independently for each channel. This avoids the difficulty of determining the subset of data on which the correction is applied. For each channel, j , bad segments are determined based on $BCT(j,:)$; data from this segments are concatenated (including all electrodes!); and PCA is applied on it. The components explaining most of the variance are subtracted to the data of channel j . Finally, corrected segments are splice into the whole data. This is done independently for each channel.

Artifacts Correction 2: *eega_tTargetPCAxEIEEG*

```
EEG ] = eega_tTargetPCAxEIEEG( EEG, nSV, vSV 'maxtime', maxtime, 'mask', mask,  
'alltime', alltime, 'allchannels', allchannels, 'savecorrected', 1 );
```

Inputs:

EEG EEGlab structure.
nSV Number of principal components to remove from the data.
vSV Amount variance to remove from the data (it removes the first n components of the data that removes a fraction of the variance up to vSV).
maxtime Maximum duration of a segment in order to apply the correction (in seconds). Default 0.200.
mask The correction is also applied on segments around bad periods (in seconds). Default 0.
alltime Logical value defining if the correction is applied on all samples which are bad in BCT (1) or only in that ones defined as bad in BCT but for which the an artifact was defined in BT (0). Default 0.
allchannels Logical value defining if the correction is applied on all channels (1) or only in channels defined as good in BC (0). Default 0.
savecorrected Logical value defining if the corrected data will be marked as good in BCT (1) or not (0). Default 1.

Outputs:

EEG EEGlab structure with the update rejection structure.

- **eega_tInterpSpatial** This function use spherical spline to interpolate electrodes not working during an entire epoch.

Artifacts Correction 3: *eega_tInterpSpatial*

```
EEG, Interp_out ] = eega_tInterpSpatial( EEG, p_int, excbad );
```

Inputs:

EEG EEGlab structure.

Outputs:

EEG EEGlab structure with the update rejection structure.

- **eega_tInterpSpatialSegment** This function use spherical spline to interpolate segments of bad data, which do not extend during the whole epoch, for each electrode.

Artifacts Correction 4: *eega_tInterpSpatialSegment*

```
EEG, Interp_out ] = eega_tInterpSpatialSegment( EEG, p_int );
```

Inputs:

EEG EEGlab structure.

Outputs:

EEG EEGlab structure with the update rejection structure.

In order define bad segments and bad channels the following function is used. After each interpolation the data marked as good in BCT is marked as good. Therefore it is necessary to define again BT and BC based on BCT. Nevertheless, the definition of BT and BC during the interpolation does not make a big difference. BC will define which channels are spatially interpolated during the whole time, thus only very bad channels should be define as bad. BT will defined on which periods the target PCA interpolation and the spatial interpolation by segments will be applied. Because each of this functions have anyway they threshold, it is better no to be strict to define BT in this stage.

- **eega_tDefBTBC** This function defines bad channels (BC) and data points (BT) based on the percentage of bad samples and bad channels respectively.

Artifacts Correction 5: *eega_tDefBTBC*

```
EEG, BTnew, BCnew ] = eega_tDefBTBC( EEG, rLimCh, rLimSmp, aLimCh, aLimSmp,
'minGoodTime', minGoodTime, 'minBadTime', minBadTime, 'maskTime', maskTime,
'badChtotdata', badChtotdata, 'keeppre', 0, 'plot', 1 );
```

Inputs:

EEG EEGlab structure.

rLimCh relative threshold for the number of bad channels, determining that a sample is bad. It is the number of IQ from the 3rd quartile.

rLimSmp relative threshold for the number of bad samples, determining that a channels is bad. It is the number of IQ from the 3rd quartile.

aLimCh absolute lower and upper limits for the maximum number of bad channels to define a bad sample (proportion). It's a 2 elements vector with numbers between 0 and 1.

aLimSmp absolute lower and upper limits for the maximum number of bad samples to define a bad channel (proportion). It's a 2 elements vector with numbers between 0 and 1.

minGoodTime Periods with non artifact shorter than minGoodTime (in seconds) are marked as having artifacts. Default 0.

minBadTime Periods with artifacts shorter than minBadTime (in seconds) are not included. Default 0.

maskTime Samples around a bad segment marked as bad in BT (in seconds). Default 0.

badChtotdata Logical value indicating if bad channels are defined based on all the recording (1) or per epoch (0). Default 0.

keeppre Logical value indicating if keeping (1) or not (0) the previous BT and BC. Default 0.

plot Logical valued indicating if producing and output plot (1) or not (0). Default 0.

Outputs:

EEG EEGlab structure with EEG.artifacts updated.

BTnew New bad segments.

BCnew New bad channels.

7 Obtain ERP

In order to obtain the ERP the data has to be locked to an event (epoched), reference (in principle to the average), and eventually z-scored and/or a baseline remove. Removing a baseline seems to make statistics worst, while normalizing by z-score using the mean and standar desviation during a long baseline seems to improve statistics.

In order to avoid distortions of the topology ERP do not correct artifacts after referencing to the average.

- **eega_epoch** This function epochs the data relative the provided event and using the provided time interval.
- **eega_tDefBEart** This function defines bad epochs based on the percentage of data rejected and interpolated.

Artifacts Correction 6: *eega_tDefBEart*

```
EEG, isbadepoch ] = eega_tDefBEart( EEG, maxBCTS, maxBTS, maxBCS, maxCCTS,  
'relative', 1, 'keeppre', 1, 'maxcicles', 4, 'plot', 1, 'savefigure', 1);
```

Inputs:

EEG EEGlab structure.

Outputs:

EEG EEGlab structure with the update rejection structure.

- **eega_tDefBEDist** This function defines bad epochs based on the percentage of data rejected and interpolated.

Artifacts Correction 7: *eega_tDefBEDist*

```
EEG, isbadepoch ] = eega_tDefBEDist( EEG, maxMeanDist, maxMaxDist, 'relative', 1,  
'keeppre', 1, 'maxcicles', 4, 'plot', 1, 'savefigure', 1 );
```

Inputs:

EEG EEGlab structure.

Outputs:

EEG EEGlab structure with the update rejection structure.

- **eega_refavg** This function subtracts the average across electrodes at each data point.
- **eega_normalization** This function normalizes the data of each electrode and epoch by z-score. The mean and standard deviation are calculated during the time interval indicated in the input.
- **pop_rmbase** This function subtracts a baseline to each electrode and epoch by removing the mean over the period indicated in the input.
- **eega_infoevents** This function recovers extra information from a text file about the events in the EEG structure. It can be applied on epoched and not epoched data.
- **eega_definefactors** This function defines relevant factors and its value for each epoch. The factors have to be fields of EEG.epoch based on the event information present in EEG.epoch.
- **eega_avgdata** This function calculates the average across epochs. The epochs are grouped and averaged together, based on the condition they belong, therefore the size of output EEG.data is electrodes x samples x conditions. Conditions are defined based on a table provided in the input, which indicates for each condition the value of the relevant factors. Consider that this function needs that EEG.F exists, thus it has to be run after *eega_definefactors*.